

INSTRUCTIONS FOR USE

EMESIS BAG DISPENSER

DESCRIPTION:

The Emesis Bag Dispenser is a wall-mounted unit made of durable, easy-to-clean materials. It is designed to hold Emesis bags, with a user-friendly mechanism to allow ease of access to emesis bags.

INDICATIONS

The Emesis Bag Dispenser is intended for use in any medical or healthcare facility where emesis bags are necessary for managing patient vomiting episodes, particularly in post-operative care, during chemotherapy treatments, or in emergency situations.

INTENDED USE:

The Emesis Bag Dispenser is intended to store and dispense emesis bags in healthcare settings such as hospitals, clinics, or surgical centers. It is designed to facilitate quick and efficient access to emesis bags by healthcare professionals to manage patient nausea and vomiting in clinical environments.

CONTRAINDICATIONS

This device is not designed, sold, or intended for use except as indicated.

KNOWLEDGE AND INTENDED SETTING

Device intended for use by trained healthcare personnel. Device use limited to surgical setting and operating room in a hospital or surgery center.

CLINICAL BENEFITS:

- The dispenser allows healthcare professionals to rapidly access emesis bags, reducing the time it takes to manage patient vomiting episodes.
- The dispenser provides easy access to emesis bags, which are stored in a sanitary manner, reducing the risk of contamination and infection and promoting a cleaner clinical environment for both patients and staff.

INTENDED USER

This device is intended for use by trained healthcare professionals only.

INTENDED CARE SETTINGS

Hospitals and Ambulatory Surgery Centers (ASC).

WARNINGS


- Do not overload the dispenser beyond its maximum capacity.
- Use only with compatible emesis bags recommended by Xodus Medical.

PRECAUTIONS

- Ensure the dispenser is installed properly to avoid accidental detachment from the wall.
- Periodically check the dispenser for signs of wear or damage and replace if necessary.

PREPARATION AND USE

1. Place screw through a washer with metal side on top.
2. Place neoprene side of washer against plastic and push screw through one of the smaller holes in the plastic dispenser from the inside.
3. Place another washer on the tip of the screw protruding from the other side with the neoprene side of washer against the plastic dispenser.
4. Place the assembly against the wall in the area to be attached to.
5. Screw into the wall.
6. Repeat with the second set of washers and screw.



DISPOSAL

Dispose of the dispenser and any related packaging materials in accordance with hospital protocol.

TECHNICAL SPECIFICATIONS

- Emesis Bag: Polyethylene
- Emesis Bag Ring: Polyethylene
- Dispenser: Polycarbonate

STORAGE, TRANSPORT, AND OPERATIONAL CONDITIONS

- The emesis bag dispenser should be stored in a clean, dry location at room temperature prior to use.
- Avoid exposure to extreme temperatures or humidity.



Non-sterile



xodusmedical.com/eifu

Consult instructions
for use or consult
electronic
instructions for use



Do Not Re-use



Not Made with Natural
Rubber Latex



Do Not Use If The
Package Is Damaged
Or Opened



XODUS MEDICAL
Making Surgery Safer™

702 Prominence Drive
New Kensington, PA 15068,
United States

xodusmedical.com
info@xodusmedical.com
Tel: 1+724-337-5500

FDA REGISTERED
ISO 13485 CERTIFIED



MADE IN THE
USA

LITIFUEB001 R6 10/1/2024