



# The Pink Hip Kit™

## Postless Positioning System



#72205289



Eliminate Post  
Related Injury



#72205682

[xodusmedical.com/pinkhipkit](http://xodusmedical.com/pinkhipkit)

- **Safer Surgery**

Go Postless! The Pink Hip Kit safely provides the counter-traction required to offset the force of hip distraction—securing the patient without leveraging the perineum against a rigid post.


- **Reduced Risk of Hospital Acquired Pressure Injuries (HAPIs)**

The Pink Hip Kit features a breathable, moisture-wicking, soft, enveloping, immersive, and pressure-distributing material that interfaces with sensitive patient areas. This unique material provides an ideal balance of skin-protection properties that reduce the risk of HAPIs.

- **Streamlined Procedures**

A postless approach requires shorter traction time, eases hip distraction, improves C-Arm visualization and provides a greater range-of-motion for more complete dynamic exams—all without changes to your surgical technique.

REPEATABLE  
**BEST**  
PRACTICE




Eliminate the need for custom, off-label positioning solutions

Setup time: less than one minute


Reduce cost and eliminate cross-contamination risks

### The Material Difference


Patients are immersed and enveloped in a unique, patented material that both prevents HAPI formation and patient sliding.




### Multiple Configurations Available:




**Stryker Table**  
40618STRBL




**Maquet Table**  
40624MAQ




**Arthrex Table**  
40622ARTBL



**Nuova BN Table**  
40625NU



**Smith & Nephew**  
72205289\*



**Hana Table**  
72205682\*

*\*Exclusively Distributed by Smith+Nephew*



**XODUS MEDICAL**  
Making Surgery Safer™

702 Prominence Drive  
New Kensington, PA 15068,  
United States

xodusmedical.com  
info@xodusmedical.com  
Tel: 1+724-337-5500

FDA REGISTERED  
ISO 13485 CERTIFIED



MADE IN THE  
**USA**



Follow Us on  
**Social Media**