

SURGICAL TRAY CORNER PROTECTOR




- Prevents holes, tears or rips in the corners of the wrapped surgical trays.
- Open cell foam creates a soft barrier to protect sterilization wrap from sharp corners.
- Eliminate unnecessary costs created by contaminated surgical trays.
(*Average savings \$40 per contaminated tray)
- No more double wrapping of trays.
- Cost effective solution

40461 – 100 sheets/case (4 squares per sheet)

40462 – 100 sheets/case (4 squares per sheet)



Easily Separable



Pittsburgh, PA U.S.A • Tel: 724-337-5500 Tel: 800-963-8776
Fax: 724-337-0555 Fax: 800-963-6553

www.xodusmedical.com • email: info@xodusmedical.com

