




# Electrosurgical

Pencils, Electrodes,  
and Accessories



A Complete Line  
of Superior Quality  
Electrosurgical  
Products

# Laparoscopic Tips



## Product Numbers

Size	33 cm.	37 cm.	44 cm.
Hand Control	#20080		#20145
Foot Control	#20127		#20111



Size	33 cm.	37 cm.	44 cm.
Hand Control	#20081		#20138
Foot Control			



Size	33 cm.	37 cm.	44 cm.
Hand Control	#20082		#20146
Foot Control	#20115		#20147



Size	33 cm.	37 cm.	44 cm.
Hand Control	#20106	#20158	#20150
Foot Control			



Size	33 cm.	37 cm.	44 cm.
Hand Control	#20135		
Foot Control			



Size	33 cm.	37 cm.	44 cm.
Hand Control	#20144		#20152
Foot Control	#20124		#20153

## Patient Safety First

### 100% Tested

- Industry leading, heat-resistant Kynar insulation
  - » Single-use eliminates reprocessing
  - » Eliminates cross-contamination
  - » Lowers total costs
  - » Prevents insulation damage, associated burn risk
- Compromised insulation may lead to stray currents and thermal necrosis of non-targeted tissue—which carries the potential for peritonitis, other serious injury or even death
- Every electrode individually tested (high-potential stress test)
- PTFE coatings deliver consistent coagulation and resist flaking at high temperatures



### PTFE Coating

Ultra durable coating resists eschar build up, cleans easily

### Kynar™ Insulation

Extreme heat resistant insulation withstands the rigors of electrosurgery

### Polyolefin Insulation















Minimize stray-sight burn risk with individually tested, full length insulation

### Stainless Steel Core




High quality, solid cores improve durability, and conductivity

Kynar is a registered trademark of Arkema, Inc.

# PTFE Coated Electrodes

	Product No: #20023	50 /Case
	Product No: #20062	50 /Case
	Product No: #20014	50 /Case
	Product No: #20015	50 /Case
	Product No: #20018	50 /Case
	Product No: #20012	50 /Case
	Product No: #20013	50 /Case
	Product No: #20057	50 /Case
	Product No: #20058	50 /Case
	Product No: #20055	50 /Case
	Product No: #20056	50 /Case
	Product No: #20010	50 /Case
	Product No: #20011	50 /Case
	Product No: #20017	50 /Case

# Tungsten Micro-Needles


	Product No: #20072	12 /Case
	Product No: #20070	12 /Case
	Product No: #20102	12 /Case

# Pencils with PTFE Coated Electrodes

High quality components and construction allow for precise switch control, safeguarding against accidental activation. Performance is exceptional at lower power settings, helping to minimize thermal necrosis risk. Patient safety is maximized through strict quality control and the use of premium materials. Optional preloaded, PTFE coated electrodes minimize eschar build-up.

#21850

Rocker Switch



#21800

Button Switch

All pencils are available with PTFE or Stainless Tips  
View a full list of pencil and electrode combinations  
on [xodusmedical.com](http://xodusmedical.com)



Holsters Included With All  
Electrosurgical Pencils

## Grounding Pads



- Round, latex-free, split design with and without a 10 ft (3m) cord
- Overlapping gel-line border limits aluminum edge exposure, reducing burn risk
- Compatible with generators featuring contact quality monitors

Product No.	Description	Case QTY
#20240	Split Grounding Pad with Cord	50
#20230	Split Grounding Pad without Cord	50

## Cautery Tip Cleaners



- Abrasive surface effectively removes tissue and debris with one pass
- Cleans uncoated (stainless steel) cautery tips
- Adhesive backing for universal placement
- X-ray detectable

Product No.	Case QTY
#30500	50

## EZ Covers



- Prevents cross-contamination
- Easy open/close draw string
- Easy cleanup of blood & other infectious materials

Product No.	Case QTY
#50300	50



**XODUS MEDICAL**  
Making Surgery Safer™

702 Prominence Drive  
New Kensington, PA 15068,  
United States

[xodusmedical.com](http://xodusmedical.com)  
[info@xodusmedical.com](mailto:info@xodusmedical.com)  
Tel: 1+724-337-5500

FDA REGISTERED  
ISO 13485 CERTIFIED



MADE IN THE  
**USA**



Follow Us on  
**Social Media**

LITBRES006 R4 4/2/2025