



The Pink Hip Kit™

Postless Positioning System

Eliminate Post
Related Injury



Pink Hip Kit:
Smith + Nephew
#72205289



Pink Hip Kit:
Hana
#72205682



xodusmedical.com/pinkhipkit

- **Safer Surgery**

Go Postless! The Pink Hip Kit safely provides the counter-traction required to offset the force of hip distraction—securing the patient without leveraging the perineum against a rigid post.

- **Reduced Risk of Hospital Acquired Pressure Injuries (HAPIs)**

The Pink Hip Kit features a breathable, moisture-wicking, soft, enveloping, immersive, and pressure-distributing material that interfaces with sensitive patient areas. This unique material provides an ideal balance of skin-protection properties that reduce the risk of HAPIs.

- **Streamlined Procedures**

A postless approach requires shorter traction time, eases hip distraction, improves C-Arm visualization, and provides a greater range-of-motion for more complete dynamic exams—all without changes to your surgical technique.



"I have fully adopted Pink Hip Kits as my preferred method for achieving postless hip distraction. Beyond the patient safety advantages of the surgical technique, I am confident that the Pink Hip Kit itself has a powerful, positive impact on my patients' outcomes."

-Mia S. Hagen, M.D.
Associate Professor, Department of Orthopaedics & Sports Medicine
Associate Program Director, Orthopaedic Surgery Residency
Team Physician, Washington Huskies

REPEATABLE BEST PRACTICE

Eliminate the need for custom, off-label positioning solutions




Setup time: less than one minute




Reduce cost and eliminate cross-contamination risks

The Material Difference



Patients are immersed and enveloped in a unique material that both prevents HAPI formation as well as patient sliding




Multiple Configurations Available:




Stryker Table
#40618STRBL




Maquet Table
#40624MAQ




Arthrex Table
#40622ARTBL



Nuova BN Table
#40625NU



Smith + Nephew
#72205289*



Hana Table
#72205682*

*Exclusively Distributed by Smith+Nephew



XODUS MEDICAL
Making Surgery Safer™

702 Prominence Drive
New Kensington, PA 15068,
United States

xodusmedical.com
info@xodusmedical.com
Tel: 1+724-337-5500

FDA REGISTERED
ISO 13485 CERTIFIED



MADE IN THE
USA



Follow Us on
Social Media