

Part of the Camel™ Line of Fluid Waste Management / Floor Safety Products

The Oasis™


PATENT PENDING

The industry's only anti-fatigue mat with fluid absorption properties, offering maximum protection to the OR staff and patient.

The OR can be a dangerous environment. Muscle fatigue can lead to diminished surgical performance, endangering staff and patient, while fluid build-up in the OR can present its own risks. The Oasis™ deters the onset of muscle fatigue as it helps to keep the OR clean, dry and safe from slips, falls and contaminants.

At the core of the Oasis is the same ultra absorbent material featured in the Camel® Mat — sandwiched between two layers of shock absorbing, closed-cell foam. A fluid-proof poly barrier prevents fluid from reaching the OR floor, while the non-slip backing keeps the mat securely in place.

Stacking up the Protection



Product # 80540

- Minimize the effects of muscle fatigue to OR staff
- Maximize performance – Muscle fatigue adversely impacts surgical performance
- Waste fluid management – lessen volume of surgical fluid spread throughout the OR
- Infection control – reduce exposure to contaminants
- Increase safety – minimize possibility of injury from slips and falls
- Save money – reduce OR turnover time



Call for a **FREE** sample

Pittsburgh, PA USA • Tel: 724-337-5500 • Fax: 724-337-0555
xodusmedical.com • info@xodusmedical.com

CE ISO 13485 CERTIFIED
FDA REGISTERED

latex free MADE IN THE USA

LITBRAM002 R2 2/14/18