


LLETZ Loops

Xodus Medical's LLETZ Loops are perfect for office and outpatient cervical conization.

The tungsten wire offers precision and durability at higher power settings, minimizing thermal damage to biopsy specimens.




Electrodes arrive in a protective plastic case to preserve wire integrity.



20201 Round - 10mm(W) X 8mm(H)



20202 Round - 10mm(W) X 10mm(H)




20203 Round - 15mm(W) X 12mm(H)



20204 Round - 15mm(W) X 15mm(H)



20205 Round - 20mm(W) X 12mm(H)



20206 Round - 20mm(W) X 15mm(H)



20210 Square - 10mm(W) X 10mm(H)



XODUS MEDICAL
Making Surgery Safer™



**MADE IN THE
USA**

Pittsburgh, PA USA • Tel: 724-337-5500 • Fax: 724-337-0555
xodusmedical.com • info@xodusmedical.com

ISO 13485 CERTIFIED
FDA REGISTERED