


Attaches directly to 5mm and 12mm trocars for quick and thorough cleaning

Scratch-free cleaning sponges

High-performance batteries last up to five hours



Superior laparoscopic cleaning restores optimum visibility

Thermostatically controlled reservoir keeps anti-fog solution just above body temperature

Dual layer sponge envelops and protects the lens for thorough, equipment-friendly cleansing



Product #30977

All Encompassing Solution

This comprehensive kit ensures simple, safe, superior lens cleaning from beginning to end of your MIS cases. The kit includes: One See Sharp unit, one white balancing microfiber cloth, two bottles of anti-fog solution, one 5mm and one 12mm safety-ring trocar swab.

- **EFFICIENCY** Quick, convenient access to laparoscope cleaning and defogging right on the surgical field, letting you focus on patient safety.
- **EFFICACY** Heated cleaning solution and scratch-free cleaning media restore optimum visibility.
- **TOTAL CONTROL** Recessed switch gives you the power to control when the unit is on and off, conserving the five-hour battery life.