

The Pink Hip Kit™

Go Postless!

■ Prevents Pressure Related Injury to the Perineum


By removing the perineal post, The Pink Hip Kit eliminates risk of nerve and tissue injury to the groin area and significantly reduces risk of HAPI formation

■ Streamlines Procedure

Our postless approach requires shorter traction time, makes distraction easier, and improves C-arm visualization to allow for a more complete dynamic exam without any change to surgical technique

■ New Standard of Care

The Pink Hip Kit affords a greater range of motion during surgery that makes complex procedures easier and significantly reduces post-operative pain and recovery times



*Exclusively distributed by Smith+Nephew

Pink Hip Kit Includes:

- 1 Pink Hip Pad
- 1 Pink Pad Lined Head Positioner
- 1 Two-Piece Body Strap
- 1 Pair of Armboard Pads
- 1 Pair of Armboard Straps
- 1 Lift Sheet



XODUS MEDICAL
Making Surgery Safer™

ISO 13485 CERTIFIED
FDA REGISTERED

702 Prominence Drive
New Kensington, PA 15068

Tel: 724-337-5500
Fax: 724-337-0555

xodusmedical.com
info@xodusmedical.com

CE

MADE IN THE
USA

latex
free

LITBRPP014 R4 3/3/2020