

# SAFETY SCALPEL

**Eliminate Transmission of Blood-Borne Pathogens**


**Comply with OSHA Requirement 1910.1030(d)(2)(i)**

“Engineering and work practice controls shall be used to eliminate or minimize employee exposure.”

Protective safety cartridges completely cover blades

Cartridge ejection button adds an extra layer of safety

Reusable handle can be cleaned and sterilized like a traditional handle



Same full blade exposure as a traditional scalpel

Safety cartridge locks into place

Stainless steel handle has same length, weight, shape, center of gravity and “feel” of a traditional scalpel

Safety Blade Cartridges



**90100**  
Stainless  
#10 Blade



**90101**  
Carbon  
#10 Blade



**90110**  
Stainless  
#11 Blade



**90111**  
Carbon  
#11 Blade




**90150**  
Stainless  
#15 Blade



**90151**  
Carbon  
#15 Blade

Handles

- 90000** - #3 Handle - Stainless Steel
- 90001** - #3 Long Handle - Stainless Steel
- 90007** - #7 Handle - Stainless Steel



**XODUS MEDICAL**  
Making Surgery Safer™

ISO 13485 CERTIFIED  
FDA REGISTERED

Pittsburgh, PA USA • Tel: 724-337-5500  
Fax: 724-337-0555

xodusmedical.com  
info@xodusmedical.com

CE  
0459

MADE IN THE USA  
latex free

LITBRIP003 R3 11/26/18

# SINGLE-USE SAFETY SCALPEL


Eliminate Transmission of Blood-Borne Pathogens

**Comply with OSHA Requirement 1910.1030(d)(2)(i)**

“Engineering and work practice controls shall be used to eliminate or minimize employee exposure.”

Protective safety cartridge completely covers the blade

Handle has same length and shape of a traditional scalpel



Same full blade exposure as a traditional scalpel

Safety cartridge locks into place

## Full Line of Single-use Safety Scalpels



90210 - Single Use Safety Scalpel #10 Stainless Blade





90211 - Single Use Safety Scalpel #11 Stainless Blade



90215 - Single Use Safety Scalpel #15 Stainless Blade

**Ideal for Emergency Rooms  
and Cath Labs**

Available Non-Sterile in your custom packs  
(i.e. C-Section Kits or Minor Procedure Kits)



 **XODUS MEDICAL**  
Making Surgery Safer™

ISO 13485 CERTIFIED  
FDA REGISTERED

Pittsburgh, PA USA • Tel: 724-337-5500  
Fax: 724-337-0555

xodusmedical.com  
info@xodusmedical.com

CE  
0459

MADE IN THE  
**USA** 

LITBRIP004 R2 2/5/18