

STRETCHER PADS

Replacement Transport Pads

- Available in multiple sizes to fit any stretcher
- Polyurethane coated nylon cover material prevents body fluid from penetrating foam core
- Guarantee secure placement with large adhesive on back
- Fire retardant



Pittsburgh, PA U.S.A • Tel: 724-337-5500
Fax: 724-337-0555

www.xodusmedical.com
info@xodusmedical.com



ISO 13485 CERTIFIED
FDA REGISTERED



LITBRSP001 R2 2/12/18