




The **Pink** Hip Kits™



**Eliminate Post Related
Injury**

xodusmedical.com/pinkhipkit

■ **Safer Surgery**

Go Postless! The Pink Hip Kit safely provides the counter-traction required to offset the force of hip distraction—securing the patient without leveraging the perineum against a rigid post.

■ **Reduced Risk of Hospital Acquired Pressure Injuries (HAPIs)**

The Pink Hip Kit features a breathable, moisture-wicking, soft, enveloping, immersive, and pressure-distributing material that interfaces with sensitive patient areas. This unique material provides an ideal balance of skin-protection properties that reduce the risk of HAPIs.

■ **Streamlined Procedures**


A postless approach requires shorter traction time, eases hip distraction, improves C-Arm visualization and provides a greater range-of-motion for more complete dynamic exams—all without changing to your surgical technique.

**REPEATABLE
BEST
PRACTICE**

Eliminate the need for custom, off-label positioning solutions.




Setup time: less than one minute.




Reduce cost and eliminate cross-contamination risks.

The Material Difference



Patients are immersed and enveloped in a unique, patented material that both prevents HAPI formation and patient sliding.





Multiple Configurations Available:





Stryker Table
40618STRBL





Smith & Nephew
72205286




Arthrex Table
40616ART



Nuova BN Table
40623NU



Maquet Table
40624MAQ



Hana Table
HIP40615H



702 Prominence Drive
New Kensington, PA 15068

Tel: 724-337-5500
Fax: 724-337-0555

xodusmedical.com
info@xodusmedical.com



FDA REGISTERED
ISO 13485 CERTIFIED



MADE IN THE
USA



LITBRPP022 R1 8/24/2021